



### 3<sup>rd</sup> Generation Lentiviral Vector System

#### “Purer” GMP Plasmid Manufacturing




- Global supply: both **CN** and **US** plasmid facility
- Licensed strain is **clear sourced without IP issue**
- Ready-to-use **LentiHelper™** plasmids with **FDA DMF Registration** (MF25552)

#### Viral Vector Packaging ——— Cell Therapy




- Licensed adherent cell line: **PowerA™-293T** with CDMO authority
- Proprietary suspension cell line: **PowerS™-293T** with **FDA DMF registration** (MF29374/0)

**Solid Track Records with over 60 Global IND Approvals**  
**Proven Solution with 500+ Clinical Manufacturing batches**




## 6-Month One-stop Solution from Plasmid to Lentiviral Vector




## “Purer” GMP Plasmid Manufactured & Delivered in Less than Half of Market Lead-time

□ High Density Fermentation: up to 1g/L



□ 95% supercoiled plasmid




□ Impurity removal, 1 EU/mg endotoxin



- Host DNA: <1% w/w VS <5% w/w
- Host Protein: <1% w/w VS 2% w/w
- Host RNA: <1% w/w VS <5% w/w
- Endotoxin: low as 1 EU/mg VS 100EU/mg mkt benchmark (customized item)

## LVV: High Titer and Yield

□ Stable and High Yield (1E7~1E8 TU/ml)



□ LVV produced by PowerS™-293T cell line shows higher titer than other cell lines



## Company Information

**Address** / ProBio Inc. Building 9, 311 Pennington Rocky Hill Rd, Pennington, NJ 08534

**Phone** / 1-732-885-9188 **Toll-Free** / 1-877-436-7274 **Fax** / 1-732-210-0262

**Contact** / [cdmo.us@probiocdm.com](mailto:cdmo.us@probiocdm.com) [cdmo.eu@probiocdm.com](mailto:cdmo.eu@probiocdm.com)  
[cdmo.apac@probiocdm.com](mailto:cdmo.apac@probiocdm.com)